

Letter

AIT test has no problem in the detection of anti-ribosomal P – authors' response

Michael Mahler¹ and Marvin J Fritzler²

¹Dr Fooke Laboratorien GmbH, Mainstraße 85, 41469 Neuss, Germany

²Faculty of Medicine, HRB410B, University of Calgary, 3330 Hospital Dr NW, Calgary, Alberta, Canada T2N 4N1

Corresponding author: Marvin J Fritzler, fritzler@ucalgary.ca

Published: 25 August 2009

This article is online at <http://arthritis-research.com/content/11/4/411>

© 2009 BioMed Central Ltd

Arthritis Research & Therapy 2009, **11**:411 (doi:10.1186/ar2776)

See related letter by Jearn and Kim, <http://arthritis-research.com/content/11/3/407>, and related research by Mahler *et al.*, <http://arthritis-research.com/content/10/6/R131>

In a recent article, Jearn and Kim assessed the ability of indirect immunofluorescence on a novel human macrophage cell line (IT-1) in an autoimmune target (AIT) assay to predict the presence of anti-ribosomal P (anti-Rib-P) antibodies [1], based on a previously described cytoplasmic and/or nucleolar staining pattern [2]. Despite the clearance of the AIT test by the South Korea Food and Drug Administration, to the best of our knowledge this test is not widely used as a screening method – neither is it approved by the Food and Drug Administration (FDA) USA or is it a Communauté Européenne (CE) certified method for the detection of anti-nuclear antibodies. By contrast, many laboratories rely on HEp-2 cell substrates that are widely available as FDA and CE approved diagnostic kits from a variety of manufacturers.

In their study, Jearn and Kim identified sera with anti-Rib-P reactivity based on a double immunodiffusion (ID) assay, and then used these sera to determine the positivity and staining patterns in the AIT test. This method raises concerns because it is generally known that immunodiffusion is a relatively insensitive technique for the detection of human autoantibodies when compared with other techniques such as ELISA and western immunoblotting (WB).

In an example relevant to the anti-Rib-P system, Bonfa and colleagues clearly demonstrated that immunodiffusion has very low sensitivity for the detection of anti-Rib-P antibodies [3]. In their study, only 14% of anti-Rib-P-positive samples identified by WB were detected by immunodiffusion. Another study showed that ID exhibits significantly lower sensitivity compared with a multiplex system: 63/130 (48.5%) samples with at least one positive extractable nuclear antibody result were negative by ID, and all 11 sera with anti-Rib-P reactivity

did not show a typical cytoplasmic staining pattern by indirect immunofluorescence [4]. While the aforementioned studies used HEp-2 cells as the substrate, there would be no reason to believe that the AIT assay would perform remarkably better.

In the context of contemporary diagnostic assays such as WB, ELISA or multiplexed addressable laser bead immunoassays, therefore, the data presented by Jearn and Kim do not support the claim that the AIT test represents a sensitive and reliable method for the detection of anti-Rib-P antibodies. A systematic study is mandatory to compare the sensitivity and specificity of the AIT test for the detection of anti-Rib-P antibodies and other antibody specificities using a variety of contemporary diagnostic platforms, not just ID. In addition to anti-Rib-P, such a study could also include the evaluation of anti-Jo-1, anti-SS-A (Ro60) and anti-Ro52 antibodies, which have been shown to produce false negative results in indirect immunofluorescence on HEp-2 cells [2,5].

Competing interests

MM is employed at Dr Fooke Laboratorien GmbH selling autoantibody ELISAs. MJF receives honoraria from ImmunoConcepts Inc. (Sacramento, CA, USA) for consulting services.

References

1. Jearn L-H, Kim T-Y: **AIT test has no problem in the detection of anti-ribosomal P.** *Arthritis Res Ther* 2009, **11**:407.
2. Mahler M, Ngo JT, Schulte-Pelkum J, Luettich T, Fritzler MJ: **Limited reliability of the indirect immunofluorescence technique for the detection of anti-Rib-P antibodies.** *Arthritis Res Ther* 2008, **10**:R131.
3. Bonfa E, Gaburo Júnior N, Tavares AV, Cossermelli W: **Comparison of five methods for the detection of antiribosomal P protein antibody.** *Braz J Med Biol Res* 1994, **27**:637-643.
4. Desplat-Jego S, Bardin N, Larida B, Sanmarco M: **Evaluation of**

AIT = autoimmune target; anti-Rib-P = anti-ribosomal P; ELISA = enzyme-linked immunosorbent assay; ID = immunodiffusion; WB = western immunoblotting.

the BioPlex 2200 ANA screen for the detection of antinuclear antibodies and comparison with conventional methods. *Ann N Y Acad Sci* 2007, **1109**:245-255.

5. Kidd K, Cusi K, Mueller R, Goodner M, Boyes B, Hoy E: **Detection and identification of significant ANAs in previously determined ANA negative samples.** *Clin Lab* 2005, **51**:517-521.